recon_PUP_FF_SOFC_hires3_S2

February 8, 2019

1 Overview

This notebook provides detailed steps for the reconstruction (including data reduction and morphology reconstruction) of proejct PUP_FF_SOFC_hires3_S2, which was taken place back in 2015.

NOTE:

* Currently, only the first layer (out of seven valid layers) is analyzed. * The choice of the first layer is based on the fact that it contains the tip, the poristy of which has significant impact on the FC performance.

1.1 Prepation

Import modules and define user define functions.

```
In [1]: # -- import modules
        import tomopy
        import gc
        import numpy
                       as np
        import scipy
                       as sp
        import numexpr as ne
        import matplotlib.pyplot as plt
        from tqdm import tnrange, tqdm_notebook
        from tqdm import trange
                               import read_aps_1id_metafile
        from dxchange.reader
        from dxchange.reader
                               import read_tiff_stack
                               import write_aps_1id_report
        from dxchange.writer
        from dxchange.writer
                               import write_tiff_stack
        from dxchange.writer
                               import write_vtr
        from tomopy.prep.alignment import find_slits_corners_aps_1id
        from tomopy.prep.alignment import detector_drift_adjust_aps_1id
        from tomopy.prep.alignment import remove_slits_aps_1id
        from tomopy.prep.normalize import minus_log
        from tomopy.misc.npmath
                                   import discrete_cdf
        from tomopy.misc.npmath
                                   import gauss1d
```

```
from scipy.signal import medfilt
        from scipy.signal import medfilt2d
        # from scipy.signal import convolve
        from scipy.optimize import curve_fit
        from scipy.ndimage import gaussian_filter
        from scipy.ndimage import shift
        from scipy.ndimage import convolve
        from scipy.special import erf, erfc
        # --
        import ipywidgets as widgets
        from ipywidgets import interact, interactive, fixed, interact_manual, IntProgress
        from IPython.display import display
        # --
        from matplotlib import animation, rc
        from IPython.display import HTML, Image
/home/chenz3/.pyenv/versions/anaconda3-5.2.0/lib/python3.6/site-packages/skimage/transform/_was
  warn("The default mode, 'constant', will be changed to 'reflect' in "
In [2]: # -- svd function
        def svd_enhance(img, eigen_cut=20):
            U, S, V = np.linalg.svd(img, full_matrices=True)
            eigen_cut = min(eigen_cut, U.shape[1], V.shape[0])
            return np.dot(U[:,:eigen_cut]*S[:eigen_cut], V[:eigen_cut,:])
In [3]: # -- quick lambda function
        val_atPercent = lambda ar, p: np.sort(ar.flatten())[int(np.prod(ar.shape)*p)]
        wgt_histequal = lambda ar: (np.sort(ar.flatten()).searchsorted(ar) + 1)/np.prod(ar.sha
        wgts_binned = lambda wgts, bins: np.int64(np.floor(wgts * bins)).astype(wgts.dtype)
```

1.2 Metadata processing

Process/Parse the metadata (auto-generated during scan) to extract

- beam status (various stability)
- frame quality
- connectivity bewteen rotation angle (ω) and the image sequence number
- group layers into dataframes

```
In [4]: # ---- CONFIG BLOCK ---- #
    # -- change the following based on the sample
    metaf1 = "PUP_FF_SOFC_hires3_S2_TomoStillScan.dat"
    sampleImgName = "PUP_FF_SOFC_hires3_S2/PUP_FF_SOFC_hires3_S2_000008.tif"
```

```
= "PUP_FF_SOFC_hires3_S2_results/PUP_FF_SOFC_hires3_S2"
        write_base
                       = "PUP_FF_SOFC_hires3_S2_report/scanSummary.pdf"
        reportfn
        # -- reconstruction config
        recon_config = {'algorithm': 'gridrec',
                         'filter'
                                     : 'hann',
                        }
        # -- build data frame and generate report
        df_layers = write_aps_1id_report(read_aps_1id_metafile(metaf1), reportfn)
        layerIDs = df_layers['layerID'].unique()
        print(f"Unique layers are: {layerIDs}")
Unique layers are: [ 6 7 8 9 10 11 12]
   manual select the first layer, which contains the tip
In [5]: layerID = 6
        df_layern = df_layers[df_layers['layerID'] == layerID]
        df_layern.describe()
Out[5]:
                        num
                                     nSeq
                                                               f0mega
                                                                              Omega
                             3631.000000
               3631.000000
                                           3.631000e+03
                                                          3631.000000
                                                                       3631.000000
        count
        mean
                1786.165244
                             1816.000000
                                           1.500000e-01
                                                            -1.487194
                                                                          -1.487225
                1048.041119
                             1048.323741
                                           1.015994e-14
                                                           104.810435
                                                                        104.810434
        std
                                           1.500000e-01
                                                          -180.000000
                   1.000000
                                1.000000
                                                                       -180.000000
        min
        25%
                 878.500000
                              908.500000
                                           1.500000e-01
                                                           -92.250000
                                                                        -92.250050
        50%
               1786.000000
                             1816.000000
                                           1.500000e-01
                                                            -1.500000
                                                                         -1.500080
        75%
               2693.500000
                             2723.500000
                                           1.500000e-01
                                                            89.250000
                                                                         89.249900
                             3631.000000
                                           1.500000e-01
                                                           180.000000
               3601.000000
                                                                        180.000000
        max
                       ICint
                                      IC-E1
                                                     IC-E2
                                                                   IC-E3
                                                                                 IC-E4
                                              3631.000000
               3.631000e+03
                               3631.000000
                                                             3631.000000
                                                                           3631.000000
        count
        mean
               6.500000e-01
                              76380.527403
                                             78185.773891
                                                            73212.141283
                                                                            321.169375
               4.441504e-14
                               4045.376572
                                                             3879.211524
                                                                             15.836482
        std
                                              4142.571960
        min
               6.500000e-01
                                  0.000000
                                                 0.000000
                                                              419.000000
                                                                            284.000000
        25%
               6.500000e-01
                              76262.000000
                                             78056.000000
                                                            72992.000000
                                                                            309.000000
        50%
               6.500000e-01
                              76515.000000
                                             78320.000000
                                                            73298.000000
                                                                            320.000000
        75%
               6.500000e-01
                              76965.000000
                                                                            334.000000
                                             78796.000000
                                                            73905.000000
        max
               6.500000e-01
                              78201.000000
                                             81453.000000
                                                            75524.000000
                                                                            363.000000
                                                              encoderX
                                                                              seqNum
                           bpEds
                                       LN2temp
                                                  SRcurrent
                     3631.000000
                                  3631.000000
                                                3631.000000
                                                                3631.0
                                                                        3631.000000
        count
                      -13.517194
                                     84.220476
                                                 102.049297
                                                                        1815.026990
                                                                   0.0
        mean
```

std	 1.058316	0.003734	0.174000	0.0	1048.324049
min	 -13.574900	84.211000	101.639000	0.0	1.000000
25%	 -13.573100	84.218000	101.927000	0.0	907.500000
50%	 -13.572800	84.221000	102.052000	0.0	1815.000000
75%	 -13.572400	84.223000	102.170000	0.0	2722.500000
max	 7.125490	84.230000	102.508000	0.0	3631.000000

	layerID	TBS	VBS	BLS	BLT
count	3631.0	3631.000000	3622.000000	3631.000000	3631.000000
mean	6.0	0.987991	inf	2.105510	0.006412
std	0.0	0.052350	NaN	0.110042	0.039903
min	6.0	0.005654	0.949812	0.000000	0.004191
25%	6.0	0.985021	0.976874	2.107535	0.004288
50%	6.0	0.989150	0.977052	2.112055	0.004315
75%	6.0	0.997342	0.977191	2.114556	0.004344
max	6.0	1.019190	inf	2.147057	0.768496

[8 rows x 23 columns]

1.3 Data reduction

The raw data is way to large for regular computer server to handle, therefore data reduction is necessary before further analysis.

first, generate necessary folders for subsequent analysis

```
In [8]: %%bash
```


```
mkdir -p PUP_FF_SOFC_hires3_S2_report
mkdir -p PUP_FF_SOFC_hires3_S2_results
mkdir -p PUP_FF_SOFC_hires3_S2_runtime
```

=> quick summary of background images for layer_6:

Pre-flat image: (10, 2048, 2048) as uint16

second, background normalization using collected dark and white field images

```
Post-flat image: (10, 2048, 2048) as uint16
        Post-dark image: (10, 2048, 2048) as uint16
In [10]: # interactively view the flat field image
         # 10 dark, 10 white_pre, 10_white_post
         def viz_flat(n):
             fig, axes = plt.subplots(2, 3, figsize=(15, 10),)
             axes[0,0].imshow(flat_pre[n,:,:], 'jet')
             axes[0,0].set_title("white field (before)")
             xx,yy = discrete_cdf(flat_pre[n,:,:].flatten(), steps=500)
             axes[1,0].plot(xx,yy)
             axes[1,0].set_yscale('log')
             axes[0,1].imshow(flat_pst[n,:,:], 'jet')
             axes[0,1].set_title("white field (after)")
             xx,yy = discrete_cdf(flat_pst[n,:,:].flatten(), steps=500)
             axes[1,1].plot(xx,yy)
             axes[1,1].set_yscale('log')
             axes[0,2].imshow(dark_pst[n,:,:], 'jet')
             axes[0,2].set_title("dark field")
             xx,yy = discrete_cdf(dark_pst[n,:,:].flatten(), steps=500)
             axes[1,2].plot(xx,yy)
             axes[1,2].set_yscale('log')
         # -- use the template
         interactive_plot = interactive(viz_flat, n=(0, flat_pre.shape[0]-1)) # still field
         output = interactive_plot.children[-1]
         interactive_plot
interactive(children=(IntSlider(value=4, description='n', max=9), Output()), _dom_classes=('wie
  Static view (for non-interacitve notebook mode)
In [11]: viz_flat(4)
```



now, it is time to read in the raw projections


NOTE: careful with the memory usage (about **30.0G** per layer).

```
interactive_plot = interactive(viz_proj, n=(0, projs.shape[0]-1)) # still field
output = interactive_plot.children[-1]
interactive_plot
```

interactive(children=(IntSlider(value=1800, description='n', max=3600), Output()), _dom_classe

static view of proejctio images (for non-interactive notebook session)

In [14]: viz_proj(1800)




since the majority of the image is empty, it is better to isolate the region with sample to reduce total data size.

After several interactive visual inspections, the following cropping limit is selected for this layer.

• Layer_6: img[750:1150, 800:1300]

1.3.1 Crop down images

```
In [20]: # remove the background using both white field and dark field
        projs_clean = np.zeros(projs.shape)
         for n in tqdm_notebook(range(projs.shape[0])):
             img = np.copy(projs[n,:,:]).astype(float)
             img = np.where(img-bg_dark>0, img-bg_dark, 0)
             bg = np.where(bg_white-bg_dark>0, bg_white-bg_dark, 1)
             img = img/bg
             img = -np.log(np.where(img>0, img, 1))
             projs_clean[n,:,:] = np.where(img>0, img, 0)
HBox(children=(IntProgress(value=0, max=3601), HTML(value='')))
In [21]: # Interactive view of the projection images
         def viz_proj(n):
             fig, axes = plt.subplots(1, 2, figsize=(10, 5))
             axes[0].imshow(projs_clean[n,:,:], 'gray')
             xx,yy = discrete_cdf(projs_clean[n,:,:].flatten(), steps=500)
             axes[1].plot(xx,yy)
         # -- use the template
         interactive_plot = interactive(viz_proj, n=(0, projs.shape[0]-1)) # still field
         output = interactive_plot.children[-1]
         interactive_plot
interactive(children=(IntSlider(value=1800, description='n', max=3600), Output()), _dom_classe
In [22]: # static view
        viz_proj(1800)
```



save the reduced projections and restart kernel (manual)

```
In [23]: np.save('PUP_FF_SOFC_hires3_S2_runtime/layer6_projs_reduced', projs_clean)
```

1.4 Tomography reconstruction

With the reduced data set, it is now possible to perform the reconstruction.

1.4.1 Counter sample jittering

Due to the experiment setup, significant sample jittering (particularly horizontal one) is present in the data set. This poses a huge issue to the tomography reconstruction process which requires a common rotation axis for all image. Consequelty, a phase correlation based correction process, termed "iterative pair-wise adjustment", is used here to counter the **horizontal** sample jittering by forcing the rotation axis of all iamge to the center of the image.

Note: Although visible vertical sample jittering is also observed, there is currently no reliable way to counter its impact due to lack of common reference feature.

```
In [7]: # load the reduced data set

projs = np.load('PUP_FF_SOFC_hires3_S2_runtime/layer6_projs_reduced.npy')

for n in tqdm_notebook(range(projs.shape[0])):
        projs[n,:,:] /= projs[n,:,:].max()

print(projs.shape)

HBox(children=(IntProgress(value=0, max=3601), HTML(value='')))
```

```
(3601, 400, 500)
In [10]: # -- center the data using iterative pairwise adjustment
             {\tt n\_iters}
                            = 10
             rot_centers = []
             for i in tqdm_notebook(range(n_iters), desc='iter'):
                   # adjust the images
                   for n_img in tqdm_notebook(range(int(dn)), desc='adjust'):
                         # locate the pair-wise center through PC
                         _cnt = tomopy.find_center_pc(projs[n_img,:,:],
                                                                    projs[n_img+dn,:,:],
                                                                    rotc_guess=projs.shape[2]/2,
                         # adjust the pair images
                         # NOTE:
                               make sure have order=1 to force using linear interpolation to avoid negat
                         val_shift = np.array([0, projs.shape[2]/2 - _cnt])
                         projs[n_img,:,:] = shift(projs[n_img,:,:],      val_shift, mode='constant',
                         projs[n_img+dn,:,:] = shift(projs[n_img+dn,:,:], val_shift, mode='constant',
                   # calculate the average rotation center
                   tmp = []
                   for n_img in tqdm_notebook(range(int(dn)), desc='calcAvg'):
                         tmp.append(tomopy.find_center_pc(projs[n_img,:,:],
                                                                          projs[n_img+dn,:,:],
                                                                          rotc_guess=projs.shape[2]/2,
                                        )
                   # get the average rotation center
                   _cnt = np.average(tmp)
                   # if the average rotation center has been adjusted to middle, stop
                   if abs(_cnt - projs.shape[2]/2)<0.1:</pre>
                         break
HBox(children=(IntProgress(value=0, description='iter', max=10, style=ProgressStyle(description)
HBox(children=(IntProgress(value=0, description='adjust', max=1800, style=ProgressStyle(description='adjust', max=1800, style=ProgressStyle(description='adjust')
HBox(children=(IntProgress(value=0, description='calcAvg', max=1800, style=ProgressStyle(descr
```

During the analysis, it was noiticed that there are three corrupted frames, which will add necessary noise to the reconstruction results. Therefore, the pairs containing corrupted frames are excluded from the reconstruction process.


```
In [11]: # -- gather the horizontal jitter profile (rotation center distribution)
        dict_rotcnts = {}
         thetas = np.radians(df_layern.loc[df_layern['type']=='still', 'Omega'].values)
        dn = int(np.pi/(thetas[1] - thetas[0])) # assuming eagual spacing along omega
        rot centers = []
         for n_img in tqdm_notebook(range(0, int(dn))):
             rot_centers.append(tomopy.find_center_pc(projs[n_img, :,:], projs[n_img+dn,:,:]
        dict_rotcnts['iterPairWiseAdjust'] = rot_centers + rot_centers
         # -- find the good/bad frames
        rot_cnts = np.array(dict_rotcnts['iterPairWiseAdjust'])
         idx_goodimg = np.where(abs(rot_cnts -250)< 1)[0]
         idx_badimg = np.where(abs(rot_cnts -250)>=1)[0]
        projs_good = projs[idx_goodimg,:,:]
        thetas good = thetas[idx gooding]
        rot_cnt = np.average(rot_cnts[idx_goodimg])
        print(projs_good.shape, thetas_good.shape, rot_cnt, np.average(rot_cnts))
HBox(children=(IntProgress(value=0, max=1800), HTML(value='')))
(3594, 400, 500) (3594,) 249.99763494713412 249.98625
```

1.4.2 counter rotation axis tilt

Prelimenary analysis (see below) indicate that the rotation axis is tilted, which is most likely caused by the mis-alignment between FF detector and the sample rotation stage. In theory, the most accurate correction would be in-plane rotation to bring the 180°-pair into a leveled position. However, due to lack of reference points, the exact rotation cannot be deterministically located. Therefore, an finite steps apprixmation is used where different rotation center is used for thin slice of data (two sinograms at a time).

```
plt.plot(cnts, label=f"{nimg}", linewidth=0.5, alpha=0.5)
plt.yscale('symlog')
plt.legend()
```

Out[12]: <matplotlib.legend.Legend at 0x7f4718c83940>



```
In [ ]: # use slab of 16 pixels
        # [32, 16, 8, 4, 2,]
        size_slab = 2
                  = np.zeros((projs.shape[1], projs.shape[2], projs.shape[2]))
        recon
        rotcnts
                  = []
        for nrow in tqdm_notebook(range(0, projs.shape[1], size_slab), desc='slab'):
            # locate the rotation center
            _projs = projs[:,nrow:nrow+size_slab,:]
            rotcnt = []
            for nimg in tqdm_notebook(range(dn), desc=f'r{nrow}_omega'):
                _rotcnt.append(tomopy.find_center_pc(_projs[nimg,
                                                     _projs[nimg+dn,:,:],
                                                     rotc_guess=projs.shape[2]/2,
                                                    ),
            _rotcnt = np.array(_rotcnt + _rotcnt) # 180 -> 360
```


```
recon[nrow:nrow+size_slab,:,:] = tomopy.recon(projs_good[:,nrow:nrow+size_slab,:],
                                                          thetas_good,
                                                          center=np.average(_rotcnt[idx_goodim;
                                                          algorithm='gridrec', filter_name='ha
            rotcnts.append(np.average(_rotcnt[idx_goodimg]))
        # removal phase artifacts
        recon = tomopy.circ_mask(recon, axis=0, ratio=0.6)
        phase_recon = tomopy.retrieve_phase(recon,
                                            pixel_size=2e-5, # 200 nm pixel size on the dete
                                            dist=10,
                                                             # the notes said 36m !?
                                            energy=50,
                                                              # 45-50 keV, according to exist
        # output to disk
        for n in tqdm_notebook(range(recon.shape[0])):
            img = np.copy(recon[n,:,:])
           pha = np.absolute(np.copy(phase_recon[n,:,:]))
           pha /= pha.max()
            img -= img.min()
            img /= img.max()
           msk = gaussian_filter(pha, sigma=20)
           recon[n,:,:] = img*msk*255
        # -- write to disk
        write_tiff_stack(recon,
                         overwrite=True,
                         fname=write_base,
                         start=0,
In [ ]: # quick sanity check
        def peekrecon(n):
            fig, axes = plt.subplots(1, 2, figsize=(20,10))
            img = np.copy(recon[n,:,:])
            axes[0].imshow(img, 'gray')
            # plot the discrete rotation center
```

perform the reconstruction

axes[1].plot(rotcnts)

```
# -- use the template
interactive_plot = interactive(peekrecon, n=(0, recon.shape[0]-1)) # still field
output = interactive_plot.children[-1]
interactive_plot
```

In [15]: # static view
 peekrecon(199)



1.5 Appendix

List versions of all the modules used in this notebook environment.

Note: Ensure reproducability of the reconstruction results.

```
def module_path(mod):
    '''Returns path to the file that module *mod* comes from.
    If it doesn't come from a file, return None.'''
    if hasattr(mod, '__file__'):
        return os.path.abspath(os.path.dirname(mod.__file__))
    else:
        return None
def from_git_repo(mod):
    '''Does the module *mod* reside in a Git repository?'''
    path = module_path(mod)
    if path:
        try:
            repo = git.Repo(path)
        except:
            return False
        else:
            return True
    else:
        return False
def git_path_sha(mod, slice=slice(0, 8, 1)):
    '''Return SHA hash for the HEAD commit for the repository
    that the module *mod* resides in.'''
    repo = git.Repo(module_path(mod))
    return repo.git_dir, repo.head.commit.hexsha[:8]
def module_version(mod):
    '''Return version string for module *mod*, or nothing if
    it doesn't have a "version" or "__version__" attribute.'''
    version = []
    if hasattr(mod, '__dict__'):
        keys = []
        for key in mod.__dict__.keys():
            if key.lower() == 'version' or key.lower() == '__version__':
                v = mod.__dict__[key]
                if isinstance(v, str):
                    version.append(v)
        if keys:
            print(mod, keys)
    if version:
        return ', '.join(version)
    else:
```

```
return ''
```

```
def find_loaded_modules(only_versioned_modules=True):
    '''Return list of loaded modules for which there is a version
    number or a Git repository commit SHA.
    Return a list of *(name, version, path_to_git_repo, git_head_sha)*,
    which has an HTML property for pretty display in IPython Notebooks.
    111
    def list_of_lists_to_HTML(lists, header_row=None):
        '''Convert a list of a list of strings to a HTML table.'''
        s = '  '
        if header_row:
           s +=  '\n\t\n\t '
            s += ''.join(['%s' % item for item in header_row])
            s += ' n t / tr 
       for inner_list in lists:
            s +=  '\n\t<\tr>\n\t\t'
            s += ''.join(['%s' % item for item in inner_list])
            s += '\n\t'
        s += '\n'
        return s
    class LoadedModules(list):
        '''Very simple wrapper for a list of lists of strings, with an attribute
        for display in IPython Notebooks.'''
        def __init__(self, *args, **kwargs):
           list.__init__(self, *args, **kwargs)
        @property
        def HTML(self):
            from IPython.display import HTML
            return HTML(
                   list_of_lists_to_HTML(
                           self, header row=['Name', 'Version', 'Path', 'SHA']))
    objs = LoadedModules()
    for i, mod in enumerate(globals().values()):
        if isinstance(mod, types.ModuleType):
            if hasattr(mod, '__name__'):
               name = mod.__name__
            else:
               name = ''
           if from_git_repo(mod):
               path, sha = git_path_sha(mod)
```

```
else:
    path = ''
    sha = ''

version = module_version(mod)

if only_versioned_modules:
    flag = version or (path and sha)

else:
    flag = True

if flag:
    objs.append([mod.__name__, version, path, sha])
    objs.sort(key=lambda r: r[0])
    return objs

In [14]: find_loaded_modules().HTML

Out[14]: <IPython.core.display.HTML object>
```